

TOTAL
● ● ● ● ●
Solution

INHEALTH
TECHNOLOGIES®



CONTENTS

Pre-Operative Education	2
"What to Expect" Brochure	2
Educational Videos	2
Surgical Solutions	2
Tracheoesophageal Puncture Set (Sterile)	2
Voice Prosthesis Placement Surgical Kits (Sterile)	2
Indwelling Voice Prosthesis (Sterile)	2
Adjustable Bi-Flanged Fistula Prosthesis (Sterile)	2
Post-Operative Care	3
Electrolarynx	3
Laryngectomy Tubes	4
Post-Operative Care Kits	5
Next Steps	7
Voice Prosthesis (Clinician Placed)	7
Special Order Voice Prosthesis (Clinician Placed)	8
Voice Prosthesis (Patient Changeable)	9
Heat & Moisture Exchange (HME)	11
Voice Prosthesis Accessories and Fitting Supplies	11
Adhesive Housings	13
Tracheostoma Buttons	14
Tracheostoma Attachments	14
And Beyond	16
Hands-Free Speech Valve	16
Skin Care	16
Stoma Protection and Accessories	17
Professional Products	19
Voice Prosthesis Puncture Kits and Management	20
RENÚ® Voice and RENÚ® Gel	23
Oto-Rhino Foreign Body Remover	25
Nasal Septal Perforation Prosthesis	26
Ordering Information	27
New Patient Enrollment Instructions	27
Free Educational Resources	27

InHealth Technologies® offers a TOTAL Solution for your total laryngectomy from pre- to post-op and beyond

TOTAL Solution serves patients at all stages of voice restoration, from products that deliver solutions for pre- and post-operative care and beyond. This continuum of care is enhanced by InHealth Technologies' customer-focused team, including leading clinical specialists, experienced patient consultants and dedicated customer service members.



● PRE-OPERATIVE EDUCATION

"What to Expect" Brochure

Available online at:
www.cutt.ly/inhealthbrochure

Educational Videos

Available at InHealth Technologies® YouTube Channel
www.cutt.ly/inhealthvideo

● SURGICAL SOLUTIONS

CAUTION: FEDERAL (US) LAW RESTRICTS THESE DEVICES TO SALE BY OR ON THE ORDER OF A PHYSICIAN

Blom-Singer® Tracheoesophageal Puncture Set (Sterile)

Designed for secondary placement of a tracheoesophageal prosthesis. Includes a disposable 16 gauge curved introducer needle with a 14 gauge Teflon® sheath, a 26" (66cm) silicone catheter with a tapered reinforced tip attached to a 24" (61cm) Teflon-coated stainless steel leader cable, and a fenestrated surgical drape with adhesive around the aperture.

Ref	Description
BE 6055	Tracheoesophageal Puncture Set (Sterile)

Blom-Singer® Voice Prosthesis Placement Surgical Kits (Sterile)

Voice prosthesis placement for primary or secondary puncture. The Blom-Singer Voice Prosthesis Placement Surgical Kits are designed to be used with Blom-Singer Classic Indwelling Voice Prosthesis, Sterile Series. Sold separately.

Ref	Description
TP 1000	16 Fr. Voice Prosthesis Placement Kit (Sterile)
TP 1001	20 Fr. Voice Prosthesis Placement Kit (Sterile)

Blom-Singer® Classic™ Indwelling Voice Prostheses (Sterile)

Designed for primary TE puncture and placement procedures. Includes insertion and cleaning accessories. Clinician placed.



Length	16 Fr.	20 Fr.	Length	16 Fr.	20 Fr.
8mm	IN 1608-SO	IN 2008-SO	12mm	IN 1612-SO	IN 2012-SO
10mm	IN 1610-SO	IN 2010-SO	14mm	IN 1614-SO	IN 2014-SO

Blom-Singer® Adjustable Bi-Flanged Fistula Prostheses (Sterile)

Designed by Eric D. Blom, Ph.D to control leakage of saliva, food, and drink around a hypopharyngeal fistula. The prosthesis is adjustable in length and comes in three flange sizes, 25mm, 38mm, and 50mm, ensuring a proper fit.

Ref	Description
FM 1003*	25mm (Small) Adjustable Bi-Flanged Fistula Prosthesis, 1 ea
FM 1004*	38mm (Medium) Adjustable Bi-Flanged Fistula Prosthesis, 1 ea
FM 1005	50mm (Large) Adjustable Bi-Flanged Fistula Prosthesis without IR Catheter
FM 1010	Blom-Singer Fistula Prosthesis Insertion and Removal Catheter

*FM 1010 Insertion and Removal Catheter supplied at no charge with FM 1003 and FM 1004.



● POST-OPERATIVE CARE

Blom-Singer® Electrolarynx EL 1000

An electrolarynx is a speaking aid that produces sound vibrations that can be conducted through the throat and mouth to speak after a laryngectomy.

InHealth Technologies® and Eric D. Blom, Ph.D are pleased to offer a high quality, cost effective electrolarynx for laryngectomees. Made in US.

Ref	Description
EL 1000	Electrolarynx Digital Speech Aid (Two 9-volt batteries included)
EL 1002	Oral Adapter
EL 1004	Battery Cap
EL 1007	Adjustable Replacement Cap Top (Only for models purchased before December 2018)



Blom-Singer® StomaSoft® Laryngectomy Tubes

Blom-Singer StomaSoft Laryngectomy Tubes are flexible, soft, and compatible with HME and other accessories from InHealth Technologies®. They feature reinforced strap holders and are available with or without fenestrations.

Size/Length	I.D./O.D.	Non-Fenestrated	Fenestrated	Rx
8/36	9.5mm/12.0mm	BE 6698	BE 6698F	■
8/55	9.5mm/12.0mm	BE 6699	BE 6699F	■
9/36	10.5mm/13.5mm	BE 6700	BE 6700F	■
9/55	10.5mm/13.5mm	BE 6701	BE 6701F	■
10/36	12.0mm/15.0mm	BE 6702	BE 6702F	■
10/55	12.0mm/15.0mm	BE 6703	BE 6703F	■
12/36	13.5mm/17.0mm	BE 6704	BE 6704F	■
12/55	13.5mm/17.0mm	BE 6705	BE 6705F	■

Blom-Singer® Laryngectomy Tubes (Sterile)

Blom-Singer Laryngectomy Tubes (Sterile) are durable, lightweight silicone tubes that are used to keep the tracheostoma open after laryngectomy. The tubes do not hold an HME or other attachments.

Size/Length	I.D./O.D.	Non-Fenestrated (Sterile)	Rx
8/36	9.5mm/12.0mm	BE 6398	■
8/55	9.5mm/12.0mm	BE 6399	■
9/36	10.5mm/13.5mm	BE 6400	■
9/55	10.5mm/13.5mm	BE 6401	■
10/36	12.0mm/15.0mm	BE 6402	■
10/55	12.0mm/15.0mm	BE 6403	■
12/36	13.5mm/17.0mm	BE 6404	■
12/55	13.5mm/17.0mm	BE 6405	■

Marpac Tracheostomy Tube Holder/Collar™

Marpac Tracheostomy Tube Holder/Collar is a single use, disposable Comfortable Collar™ used to secure the tracheostomy tube. Latex free. Adult—one size fits all.

Ref	Description
204	Tracheostomy Tube Holder/Collar

Kapitex® Trachi-Swab® Hi-Vis

The Kapitex Trachi-Swab Hi-Vis is indicated for use in patients who have a tracheostoma and use a device to maintain stoma patency. It is intended for cleaning tracheostomy tubes and inner cannulas after removal from the patient. Available in 50 packs.

Ref	Description
TRSWB3050	Trachi-Swab Hi-Vis



Blom-Singer® Post-Operative Care Kits (For Total Laryngectomy)

CAUTION: FEDERAL (US) LAW RESTRICTS THESE DEVICES TO SALE BY OR ON THE ORDER OF A PHYSICIAN

The Blom-Singer® Post-Operative Care Kits were designed with valued input from clinicians and patients. These kits contain supplies to assist patients in their transition from the hospital to home after a total laryngectomy. They are available in several configurations to meet select needs. All configurations include your choice of StomaSoft Laryngectomy Tube and HMEs as listed below.

The Post-Operative Care Kit is housed in a hard, clear storage case, which also includes information on how to obtain our Next Steps follow-up kit. The Next Steps kit includes additional sample product and educational items.

All kit configurations contain the following patient care supplies:

- 1 each 204 Marpac Tracheostomy Tube Holder/Collar™
- 1 each BE 6048 Blom-Singer® Shower Guard
- 1 each Drawstring Bag
- 1 each Mirror with Stand
- 1 each Halcyon LED Flashlight

REF	Description
PK 2100	
BE 6701	StomaSoft Laryngectomy Tube 9/55 Non-fenestrated, 1 ea
BE 1055EZ	HME Cartridge ClassicFlow Heat and Moisture Exchange, pack of 30
PK 2101	
BE 6703	StomaSoft Laryngectomy Tube 10/55 Non-fenestrated, 1 ea
BE 1055EZ	HME Cartridge ClassicFlow Heat and Moisture Exchange, pack of 30
PK 2102	
BE 6701	StomaSoft Laryngectomy Tube 9/55 Non-fenestrated, 1 ea
BE 1055EF	HME Cartridge EasyFlow Heat and Moisture Exchange, pack of 30
PK 2103	
BE 6703	StomaSoft Laryngectomy Tube 10/55 Non-fenestrated, 1 ea
BE 1055EF	HME Cartridge EasyFlow Heat and Moisture Exchange, pack of 30
PK 2106	
BE 6705	StomaSoft Laryngectomy Tube 12/55 Non-fenestrated, 1 ea
BE 1055EZ	HME Cartridge ClassicFlow Heat and Moisture Exchange, pack of 30
PK 2107	
BE 6705	StomaSoft Laryngectomy Tube 12/55 Non-fenestrated, 1 ea
BE 1055EF	HME Cartridge EasyFlow Heat and Moisture Exchange, pack of 30
PK 2104	
EL 1000	Electrolarynx, 1 ea
BE 6703	StomaSoft Laryngectomy Tube 10/55 Non-fenestrated, 1 ea
BE 1055EZ	HME Cartridge ClassicFlow Heat and Moisture Exchange, pack of 30
PK 2105	
IN 1608-SO	16 Fr. 8 mm Classic Indwelling VP Sterile Option, 1 ea
IN 1610-SO	16 Fr. 10 mm Classic Indwelling VP Sterile Option, 1 ea
IN 4069	16 Fr. Plug Insert, 1 ea
BE 6703	StomaSoft Laryngectomy Tube 10/55 Non-fenestrated, 1 ea
BE 1055EZ	HME Cartridge ClassicFlow Heat and Moisture Exchange, pack of 30



PK 2105



PK 2101



● NEXT STEPS

VOICE PROSTHESES (CLINICIAN PLACED)

Blom-Singer® Classic™ Indwelling Voice Prosthesis (Sterile)

Includes sterile prostheses for primary TE puncture and placement procedures. Includes insertion and cleaning accessories. Clinician placed.

Length	16 Fr.	20 Fr.	Rx	Length	16 Fr.	20 Fr.	Rx
8mm	IN 1608-SO	IN 2008-SO	■	12mm	IN 1612-SO	IN 2012-SO	■
10mm	IN 1610-SO	IN 2010-SO	■	14mm	IN 1614-SO	IN 2014-SO	■

Blom-Singer® Classic™ Indwelling Voice Prosthesis

Includes nonsterile voice prosthesis for secondary or replacement TEP procedures, with insertion/cleaning accessories. Clinician placed.

Length	16 Fr.	20 Fr.	Rx	Length	16 Fr.	20 Fr.	Rx
4mm	IN 1604-NS	IN 2004-NS	■	14mm	IN 1614-NS	IN 2014-NS	■
6mm	IN 1606-NS	IN 2006-NS	■	16mm	IN 1616-NS	IN 2016-NS	■
8mm	IN 1608-NS	IN 2008-NS	■	18mm	IN 1618-NS	IN 2018-NS	■
10mm	IN 1610-NS	IN 2010-NS	■	20mm	IN 1620-NS	IN 2020-NS	■
12mm	IN 1612-NS	IN 2012-NS	■				

Blom-Singer® Dual Valve™ Indwelling Voice Prosthesis

When the primary valve fails, the secondary speech valve is intended to remain viable to help prevent leakage of fluids. The Blom-Singer® Dual Valve Voice Prosthesis is available with and without large esophageal and tracheal flanges. Includes nonsterile voice prosthesis for secondary or replacement TEP procedures, with insertion/cleaning accessories.

Length	20 Fr.	20 Fr. Large Flange	Rx	Length	20 Fr.	20 Fr. Large Flange	Rx
6mm	DV 2006-NS	DV 2006-LF	■	12mm	DV 2012-NS	DV 2012-LF	■
8mm	DV 2008-NS	DV 2008-LF	■	14mm	DV 2014-NS	DV 2014-LF	■
10mm	DV 2010-NS	DV 2010-LF	■				

Blom-Singer® Advantage® Indwelling Voice Prosthesis

The Blom-Singer Advantage Indwelling Voice Prosthesis is embedded with silver oxide as a material preservative.* Includes nonsterile voice prosthesis for secondary or replacement TEP procedures, with insertion/cleaning accessories.

HARD VALVE ASSEMBLY

Length	20 Fr.	Rx
4mm	AD 2004-NS	■
6mm	AD 2006-NS	■
8mm	AD 2008-NS	■
10mm	AD 2010-NS	■
12mm	AD 2012-NS	■
14mm	AD 2014-NS	■

SOFT VALVE ASSEMBLY

Length	16 Fr.	20 Fr.	Rx
4mm	AS 1604-NS	AS 2004-NS	■
6mm	AS 1606-NS	AS 2006-NS	■
8mm	AS 1608-NS	AS 2008-NS	■
10mm	AS 1610-NS	AS 2010-NS	■
12mm	AS 1612-NS	AS 2012-NS	■
14mm	AS 1614-NS	AS 2014-NS	■

*Silver oxide, a known anti-microbial substance, has been incorporated into the silicone flap valve of the Blom-Singer Advantage Indwelling Voice Prosthesis as a material preservative.



Classic (Sterile)



Classic



Dual Valve



Dual Valve Large Flange



Hard Valve Assembly



Soft Valve Assembly

SPECIAL ORDER VOICE PROSTHESES (CLINICIAN PLACED)

—DEVELOPED BY ERIC D. BLOM, PH.D

Unique problems call for unique solutions. InHealth Technologies® offers several voice prostheses to address those needs. Clinician placed.

Special Length

Available in one millimeter incremental lengths for a more precise TE puncture fit.

	5mm	7mm	9mm	Rx
16 Fr.	IN 1605-SL	IN 1607-SL	IN 1609-SL	■
20 Fr.	IN 2005-SL	IN 2007-SL	IN 2009-SL	■

Increased Resistance

Designed for patients who experience intermittent valve opening with fluctuations in esophageal pressures, resulting in air ingestion and/or leakage.

	6mm	8mm	10mm	12mm	Rx
16 Fr.	IN 1606-IR	IN 1608-IR	IN 1610-IR	IN 1612-IR	■
20 Fr.	IN 2006-IR	IN 2008-IR	IN 2010-IR	IN 2012-IR	■

Large Esophageal Flange

Provides a seal to help manage leakage around a voice prosthesis. Also, it enables downsizing from a wider diameter prosthesis to a smaller diameter. Effectiveness depends on an optimal prosthesis length fit that achieves a snug flange-to-flange seal against the wall surface(s).

	4mm	6mm	8mm	10mm	12mm	14mm	Rx
16 Fr.	IN 1604-LEF	IN 1606-LEF	IN 1608-LEF	IN 1610-LEF	IN 1612-LEF	IN 1614-LEF	■
20 Fr.	IN 2004-LEF	IN 2006-LEF	IN 2008-LEF	IN 2010-LEF	IN 2012-LEF	IN 2014-LEF	■

LARGE ESOPHAGEAL FLANGE (SPECIAL LENGTH)

	5mm	7mm	Rx
16 Fr.	IN 1605-LESL	IN 1607-LESL	■
20 Fr.	IN 2005-LESL	IN 2007-LESL	■







Large Esophageal and Tracheal Flanges

Provides a double seal in patients with enlarged puncture tracts and tear-drop or irregularly shaped punctures. Effectiveness depends on an optimal prosthesis length fit that achieves a snug flange-to-flange seal against the wall surface(s).

	4mm	6mm	8mm	10mm	12mm	14mm	Rx
16 Fr.	IN 1604-LF	IN 1606-LF	IN 1608-LF	IN 1610-LF	IN 1612-LF	IN 1614-LF	■
20 Fr.	IN 2004-LF	IN 2006-LF	IN 2008-LF	IN 2010-LF	IN 2012-LF	IN 2014-LF	■

LARGE ESOPHAGEAL FLANGE (INCREASED RESISTANCE)


	4mm	6mm	8mm	Rx
16 Fr.	IN 1604-LEIR	IN 1606-LEIR	IN 1608-LEIR	■
20 Fr.	IN 2004-LEIR	IN 2006-LEIR	IN 2008-LEIR	■



TEP Occluders

The Blom-Singer® TEP (Tracheoesophageal Puncture) Occluder is a non-valved silicone device designed for placement in, and maintenance of, the TEP following total laryngectomy. Not a voice prosthesis.

	4mm	6mm	8mm	10mm	12mm	14mm	Rx
16 Fr.	IN 1604-LETO	IN 1606-LETO	IN 1608-LETO	IN 1610-LETO	IN 1612-LETO	IN 1614-LETO	■
20 Fr.	IN 2004-LETO	IN 2006-LETO	IN 2008-LETO	IN 2010-LETO	IN 2012-LETO	IN 2014-LETO	■



VOICE PROSTHESES (PATIENT CHANGEABLE)

Blom-Singer® Duckbill Voice Prostheses

The Blom-Singer® Duckbill is a patient changeable, one-way slit valve prosthesis.

Length	16 Fr.	Rx	Length	16 Fr.	Rx
6mm	DB 16-006	■	12mm	DB 16-012	■
8mm	DB 16-008	■	14mm	BE 6001	■
10mm	BE 6000	■	18mm	BE 6002	■



Blom-Singer® Low Pressure Voice Prostheses

The Blom-Singer® Low Pressure is a patient changeable, one-way valve prosthesis.

Length	16 Fr.	20 Fr.	Rx	Length	16 Fr.	20 Fr.	Rx
6mm	BE 6009	BE 2014	■	18mm	BE 6012	BE 2026	■
8mm	LP 16-008	LP 20-008	■	22mm	BE 6013	BE 2030	■
10mm	BE 6010	BE 2018	■	25mm	BE 6014	BE 2033	■
12mm	LP 16-012	LP 20-012	■	28mm	BE 6015	—	■
14mm	BE 6011	BE 2022	■				



Blom-Singer® Special Order Low Pressure Voice Prostheses (Patient Changeable)*

Length	16 Fr. Special Increased Resistance	16 Fr. Special Length	20 Fr. Special Length	20 Fr. Special Increased Resistance	Rx
4mm	—	SPECIAL 11-902-06	SPECIAL 11-251-06	—	■
6mm	LP 16-006-IR	—	—	LP 20-006-IR	■
8mm	LP 16-008-IR	—	—	LP 20-008-IR	■
10mm	LP 16-010-IR	—	—	LP 20-010-IR	■
12mm	LP 16-012-IR	—	—	LP 20-012-IR	■
14mm	LP 16-014-IR	—	—	LP 20-014-IR	■
16mm	—	SPECIAL 11-902-04	SPECIAL 11-251-04	—	■
18mm	LP 16-018-IR	—	—	—	■
20mm	—	SPECIAL 11-902-05	SPECIAL 11-251-05	—	■
22mm	LP 16-022-IR	—	—	—	■
25mm	LP 16-025-IR	—	—	—	■
28mm	LP 16-028-IR	—	—	—	■

*May require additional time for delivery.

Blom-Singer® HME All Day/All Night Advantages around the clock



ClassicFlow® HME



EasyFlow® HME

- 24-hour use
- MucusShield™ technology helps protect foam filter
- EasyTouch™ speech button—push to speak
- Compatible with a wide variety of accessories

HEAT & MOISTURE EXCHANGE (HME)

Research shows that long-term HME use results in less mucus¹, less coughing², and decreased risk of shortness of breath, fatigue, sleep problems, anxiety, and depression^{1,2}.

Blom-Singer® EasyFlow® HME

The foam filter of Blom-Singer EasyFlow allows you to breathe more freely to match your activity level and pulmonary needs. The EasyFlow has a non-reflective beige tint.

Ref	Description
BE 1055EF	EasyFlow HME Cartridge, pack of 30 (24-hour use)

Blom-Singer® ClassicFlow® HME

The Blom-Singer ClassicFlow HME is designed for increased moisture retention and airflow resistance.

Ref	Description
BE 1055EZ	ClassicFlow HME Cartridge, pack of 30 (24-hour use)

Blom-Singer® HumidiFilter® Holder

The Blom-Singer HumidiFilter Holder is a reusable, lightweight plastic holder that holds a HumidiFilter foam filter (assembly required).

Ref	Description
BE 1060	HumidiFilter Holder
BE 1080	HumidiFilter Foam Filters, pack of 30 (Not for 24-hour use)



VOICE PROSTHESIS ACCESSORIES AND FITTING SUPPLIES

Blom-Singer® Voice Prosthesis Sizers

The Blom-Singer Voice Prosthesis Sizer is a device used by clinicians to measure the length of the tracheoesophageal puncture to assist in selection of the appropriate voice prosthesis.

Ref	Description	Rx
BE 6036	Voice Prosthesis Sizer with 3 inserters (for 16 Fr. voice prostheses)	■
BE 2036	Voice Prosthesis Sizer with 3 inserters (for 20 Fr. voice prostheses)	■
BE 6096	Replacement Inserters (for use with both BE 6036 and BE 2036 Sizers), 10 ea	■

Blom-Singer® Dilator-Sizers

The Blom-Singer Dilator-Sizer combines two tasks into one innovative product: sizing and dilation. In one step, the clinician can dilate and measure the tracheoesophageal puncture during voice prosthesis placement procedures.

Ref	Description	Rx
BE 2051	18 Fr. Dilator-Sizer (for 16 Fr. prostheses)	■
BE 2052	22 Fr. Dilator-Sizer (for 20 Fr. prostheses)	■



1. Ackerstaff AH, Hilgers FJ, et al. Improvements in respiratory and psychosocial functioning following total laryngectomy by the use of a heat and moisture exchanger. *Annals of Otolaryngology, Rhinology & Laryngology*. 1993;102(11):878-883.
 2. Ackerstaff AH, Hilgers FJ, et al. Heat and moisture exchangers as a treatment option in the post-operative rehabilitation of laryngectomized patients. *Clinical Otolaryngology*. 1995;20:504-509.

Blom-Singer® Replacement Gel Caps

Use of the Blom-Singer Replacement Gel Caps provides a smooth, rounded shape to the esophageal flange end of the voice prosthesis, enabling easy insertion into the tracheoesophageal puncture.

Ref	Description	Rx
BE 3190	Replacement Gel Caps, 16 Fr., 90 ea	■
BE 3195*	Replacement Gel Caps, 18 Fr., 90 ea	■
BE 3290	Replacement Gel Caps, 20 Fr., 90 ea	■
BE 3295	Replacement Gel Caps, 22 Fr., 60 ea	■

*BE 3195 is not compatible with LP 3100.

Blom-Singer® Gel Cap Insertion System

The Blom-Singer Gel Cap Insertion System provides an alternative method to insert a Blom-Singer Low Pressure voice prosthesis into a gel cap. Each package contains a loading tool, pushrods and gel caps. Replacement gel caps available.

Ref	Description	Rx
LP 3100*	16/20 Fr. Low Pressure Gel Cap Insertion System	■

Blom-Singer® Flushing Device and Cleaning Brushes

For routine cleaning of a voice prosthesis as recommended by your clinician.

Ref	Description	Rx
IN 4050	Flushing Device, pack of 3	
IN 4064*	Cleaning Brush, 16/20 Fr., Small, 4mm–8mm, 3 ea	■
IN 4065*	Cleaning Brush, 16/20 Fr., Medium, 9mm–16mm, 3 ea	■
IN 4066*	Cleaning Brush, 16/20 Fr., Large, 18mm–22mm, 3 ea	■

*Fits indwelling style voice prostheses only.

Blom-Singer® Plug Inserts

The Blom-Singer Plug Insert is designed for temporary management of leakage through a voice prosthesis. Not compatible with patient changeable voice prostheses.

Ref	Description	Rx
IN 4069	16 Fr. Plug Insert (Advantage and Classic, 6mm to 20mm)	■
IN 4070	20 Fr. Plug Insert (Advantage and Classic, 6mm to 20mm)	■



ADHESIVE HOUSINGS

InHealth Technologies® presents the FitSeries—adhesive housings designed to meet your skin and stoma needs after a total laryngectomy. Available in two styles, three shapes each. Compatible with HME cartridges and ATSV II.

Blom-Singer® AccuFit® Adhesive Housings

The Blom-Singer AccuFit Adhesive Housing is designed to be our strongest adhesive option while still maintaining flexibility. It has increased elasticity and is less prone to tearing. Available in three shapes: round, oval, and oval extra.

Ref	Description
BE 6082	AccuFit Adhesive Housing, Round, pack of 30
BE 6083	AccuFit Adhesive Housing, Oval, pack of 30
BE 6084	AccuFit Adhesive Housing, Oval Extra, pack of 30

Blom-Singer® HydroFit® Adhesive Housings

The Blom-Singer HydroFit Adhesive Housing is a hydrocolloid adhesive baseplate made with a protective gel-like layer and a moisture-absorbent adhesive. The HydroFit is gentle on the skin and flexible, but maintains strong adhesion. Available in three shapes: round, oval, and oval extra.

Ref	Description
BE 6085	Hydrofit Adhesive Housing, Round, pack of 30
BE 6086	Hydrofit Adhesive Housing, Oval, pack of 30
BE 6087	Hydrofit Adhesive Housing, Oval Extra, pack of 30

Blom-Singer® TruSeal® Contour™ Low Profile Adhesive Housings

The Blom-Singer TruSeal Contour Low Profile Adhesive Housing is a lightweight, flexible material made to fit various peristomal anatomies. It has a lower profile from the standard housing, reducing its height and overall profile. Available in three shapes: round, standard, and oval.

Ref	Description
BE 6074	TruSeal Contour Low Profile Adhesive Housing, Round, pack of 30
BE 6075	TruSeal Contour Low Profile Adhesive Housing, Oval, pack of 30
BE 6076	TruSeal Contour Low Profile Adhesive Housing, Standard, pack of 30



TRACHEOSTOMA BUTTONS

Barton-Mayo™ Tracheostoma Buttons*

Developed by Medical Innovations International, Inc. in association with the Mayo Clinic to provide tracheostoma valve retention for near-total laryngectomy and total laryngectomy patients.

Ref	Size	O.D.	Length	Rx
BMB091	9	13.5mm	15mm, short	■
BMB092	9	13.5mm	22mm, regular	■
BMB093	9	13.5mm	30mm, long	■
BMB101	10	15mm	15mm, short	■
BMB102	10	15mm	22mm, regular	■
BMB103	10	15mm	30mm, long	■
BMB121	12	17mm	15mm, short	■
BMB122	12	17mm	22mm, regular	■
BMB123	12	17mm	30mm, long	■
BMB141	14	19mm	15mm, short	■
BMB142	14	19mm	22mm, regular	■
BMB143	14	19mm	30mm, long	■
BMBKIT Sizing Kit (all 12 sizes)				■
BMBKIT-SHORT Patient Sizing Kit-Short (4 short sizes only)				■

*Only available for sale in the US.



TRACHEOSTOMA ATTACHMENTS

Blom-Singer® Tracheostoma Valve Housings

The Tracheostoma Valve Housing is a reusable, flexible housing that attaches the ATSV II, HME cartridge or Shower Guard over the stoma. It is applied with separate adhesive tape discs and silicone adhesive. It is available in two sizes, standard and large. It is also offered in two different materials, standard and PVC.

Ref	Description
BE 6038	Tracheostoma Valve Housing, standard, 1 ea
BE 6039	Tracheostoma Valve Housing, large, 1 ea
BE 6038-P*	Tracheostoma Valve Housing, standard, PVC, 1 ea
BE 6039-P*	Tracheostoma Valve Housing, large, PVC, 1 ea

*May not be available in all international markets.

Blom-Singer® Adhesive Discs

A variety of double-sided tape discs are available to provide attachment of the tracheostoma valve housing.

Ref	Description	Ref	Description
BE 6041	Regular Tape Discs, standard, pack of 30	BE 6034	Heavy Duty Tape Discs, standard, pack of 30
BE 6042	Regular Tape Discs, large, pack of 30	BE 6035	Heavy Duty Tape Discs, large, pack of 30
BE 6043	Thin Tape Discs, standard, pack of 30	BE 6045	Regular Foam Discs, standard, pack of 30
BE 6044	Thin Tape Discs, large, pack of 30	BE 6046	Regular Foam Discs, large, pack of 30



● AND BEYOND

HANDS-FREE SPEECH VALVE

The Blom-Singer® Adjustable Tracheostoma Valve (ATSV II), commonly referred to as our “hands-free” valve, eliminates the need to cover the stoma with a thumb or finger to produce voice with a voice prosthesis. An optional Blom-Singer HumidiFilter® Cap may be worn with the hands-free valve to provide the benefits of a heat moisture exchange system.

Blom-Singer® ATSV II HumidiFilter® Cap and 7 Foam Filters

The ATSV II HumidiFilter Cap is a twist-on plastic cap that fits over the ATSV II. Comes with seven HumidiFilter Foam Filters. Additional replacement filters sold separately. (For use with ATSV II only.)

Ref	Description
BE 1010	ATSV II HumidiFilter Cap and 7 Foam Filters, 1 kit
BE 1030	ATSV II HumidiFilter Foam Filters, pack of 30

Blom-Singer® Adjustable Tracheostoma Valve II (ATSV II)

The Adjustable Tracheostoma Valve II (ATSV II) consists of two pieces: a Diaphragm/Faceplate and Body. The assembled two-piece valve inserts into the tracheostoma adhesive housing that is applied to the skin surrounding the stoma. A replacement Diaphragm /Faceplate is available separately.

Ref	Description	Rx
BE 8025H	ATSV II Valve (Body & Diaphragm/Faceplate), 1 kit	■
BE 8026H	ATSV II Replacement Diaphragm/Faceplate, 1 ea	■

SKIN CARE

3M™ Cavilon™ No Sting Barrier Film

The gentler way to protect skin from body fluids, adhesives and friction. Alcohol-free barrier film applies without stinging. Forms a breathable, transparent film for long-lasting protection.

Ref	Description
3344*	No Sting Barrier Film, box of 25

*May not be available in all international markets.

Skin-Prep™ Liquid Film-Forming Protective Barrier Wipes

Skin-Prep is a liquid film-forming dressing that, upon application to intact skin, forms a protective film to help reduce friction during removal of tapes and films.

Ref	Description
59420425	Skin-Prep Liquid Film-Forming Protective Barrier Wipes, box of 50



Skin Tac™ Adhesive Barrier Wipes

Skin Tac liquid adhesive barrier wipes are Latex-free, hypo-allergenic, clear, non-rubber adhesive barrier wipe.

Ref	Description
MS407W*	Adhesive Barrier Wipes, box of 50

*May not be available in all international markets.

NFS Adhesive (Silicone Adhesive)

PSA Pressure Sensitive Adhesive, silicone based adhesive for use with all silicone prosthetic devices.

Ref	Description
B-530	NFS Adhesive, 1 oz

Remove™ Adhesive Remover

Gently cleans hydrocolloid, acrylic and rubber-based adhesive residues from the skin without irritation, discomfort or trauma.

Ref	Description
59403125	Adhesive Remover, box of 50

STOMA PROTECTION AND ACCESSORIES

Cloth Stoma Protectors

A discreet and comfortable way to cover the stoma. Looks like a shirt and slips easily under other layers of clothing. Made of a soft 50/50 blend of cotton/rayon, machine washable, with an adjustable strap that fits most necks. Lengths from 47 to 53cm. Each color sold separately.

Ref/Color	Quantity
BE 6201/white	1 ea
BE 6251/black	
BE 7211/navy	

ADDvox® Stoma Cover Mock Turtlenecks (Only available in the US)

The ADDvox Stoma Cover Mock Turtleneck has inner layers of material which helps keep debris from entering the stoma. They slip under a shirt like a fashion dickey and have an outer layer of soft 100% cotton knit. Each cover is 9-1/2" long x 8-3/4" wide. Velcro fastener adjusts to fit necks of 12-1/2" to 22". Each color sold separately.

Ref/Color	Quantity
BE 6227-R2/white	1 ea
BE 6228-R2/black	
BE 6229-R2/navy	

ADDvox® Stoma Cover Scarves (Only available in the US)

The ADDvox Stoma Cover Scarf has inside layers to help filter out dust, dirt and pollen, is completely washable, and needs little care. This attractive neckline scarf looks just like a regular scarf or ascot and is available in a wide range of colors for men or women. All scarves have Velcro closures and adjust to fit necks from 17-3/4" to 22" (can be made smaller with tucks or stitching). Each color sold separately.

Ref/Color	Quantity
BE 6232-R2/black	1 ea
BE 6233-R2/white	
BE 6234-R2/navy	

ADDvox® is a registered trademark of Bruce Medical Supply.



ADDvox® Cloth Bib-Style Stoma Covers (Only available in the US)

The ADDvox soft, long-lasting bib protects the tracheostoma and lets you breathe more easily. This "bib-style" Stoma Cover will look crisp and clean after repeated laundering. Five soft layers keep dust and dirt from reaching the trachea. Fits necks up to 22". Velcro closure. Each color sold separately.

Ref/Color	Quantity
BE 6224-R2/white	1 ea
BE 6225-R2/beige	
BE 6226-R2/navy	

ADDvox® Silicone Shower Shield (Only available in the US)

The one-piece ADDvox Silicone Shower Shield will keep your airway free from water. Extremely soft and pliable, it creates a soft, comfortable and conforming seal along the irregular contours of your neck, shielding your stoma, while vents on the underside allow airflow. Made of medical grade silicone, this latex-free Silicone Shower Shield is comfortable, durable, and not affected by temperature, ozone, sunlight, soap or detergents.

Ref	Description
BE 6079-R2	Silicone Shower Shield, 1 ea

Blom-Singer® Shower Guard

The Blom-Singer Shower Guard is a one-piece, right-angle plastic shower guard that comes with a Tracheostoma Standard Adhesive Housing and 30 double-sided Adhesive Tape Discs. Also compatible with the Blom-Singer StomaSoft® Laryngectomy Tube.

Ref	Description	Rx
BE 6048	Shower Guard, 1 ea	■

Blom-Singer® Foam Stoma Protector

A simple and temporary (disposable) solution to help prevent debris from entering the stoma. Adhesive strip on one edge allows for easy attachment, access to the stoma to occlude for speech with a voice prosthesis, or for coughing. Measurement listed is approximate.

Ref	Description
BE 6220	Foam Stoma Protector, Tan, Large 2.7" L x 2.7" W (7cm x 7cm), pack of 30

ADDvox® Stoma Filters with Microporous Adhesive

ADDvox Stoma Filters offer protection against dust, dirt and pollen. Designed for single use with a hypoallergenic, non-irritating, self-adhesive microporous strip that is comfortable and easy to apply. Not recommended for cannula users. 2-1/4" wide x 2-1/2" high x 1/8" thick. Available in off-white and skin tone. Made in US.

Ref	Description
BE 6222-R2	Stoma Filters, Skin tone, pack of 60
BE 6223-R2	Stoma Filters, Off-white, pack of 60



Shower Guard Housing Tape Discs



PROFESSIONAL PRODUCTS

CAUTION: FEDERAL (US) LAW RESTRICTS THESE DEVICES TO SALE BY OR ON THE ORDER OF A PHYSICIAN

VOICE PROSTHESIS PUNCTURE KITS AND MANAGEMENT

Blom-Singer® Tracheoesophageal Puncture Set (Sterile)

Ref	Description
BE 6055	Tracheoesophageal Puncture Set

Blom-Singer® Voice Prosthesis Placement Surgical Kits (Sterile)

Ref	Description
TP 1000	16 Fr. Voice Prosthesis Placement Surgical Kit (VPP)
TP 1001	20 Fr. Voice Prosthesis Placement Surgical Kit (VPP)

Blom-Singer® Voice Prosthesis Sizers

Ref	Description
BE 6036	Voice Prosthesis Sizer with 3 inserters (for 16 Fr. voice prostheses)
BE 2036	Voice Prosthesis Sizer with 3 inserters (for 20 Fr. voice prostheses)
BE 6096	Replacement Inserters (for use with both BE 6036 and BE 2036 Sizers), 10 ea

Blom-Singer® Dilator-Sizers

Ref	Description
BE 2051	18 Fr. Dilator-Sizer (for 16 Fr. prostheses)
BE 2052	22 Fr. Dilator-Sizer (for 20 Fr. prostheses)


Blom-Singer® Replacement Gel Caps

Ref	Description
BE 3190	Replacement Gel Caps, 16 Fr., 90 ea
BE 3195*	Replacement Gel Caps, 18 Fr., 90 ea
BE 3290	Replacement Gel Caps, 20 Fr., 90 ea
BE 3295	Replacement Gel Caps, 22 Fr., 60 ea

*BE 3195 is not compatible with LP 3100.

Blom-Singer® Tracheoesophageal Puncture Dilators

Ref	Description
BE 6050	18 Fr. Tracheoesophageal Puncture Dilator (for 16 Fr. prostheses)
BE 2050	22 Fr. Tracheoesophageal Puncture Dilator (for 20 Fr. prostheses)



Blom-Singer® Post-Operative Care Kits (For Total Laryngectomy)

All kit configurations contain the following patient care supplies:

- 1 each 204 Marpac Tracheostomy Tube Holder/Collar™
- 1 each BE 6048 Blom-Singer® Shower Guard
- 1 each Drawstring Bag
- 1 each Mirror with Stand
- 1 each Halcyon LED Flashlight

NEW

Ref	Description
PK 2100	
BE 6701	StomaSoft Laryngectomy Tube 9/55 Non-fenestrated, 1 ea
BE 1055EZ	HME Cartridge ClassicFlow Heat and Moisture Exchange, pack of 30
PK 2101	
BE 6703	StomaSoft Laryngectomy Tube 10/55 Non-fenestrated, 1 ea
BE 1055EZ	HME Cartridge ClassicFlow Heat and Moisture Exchange, pack of 30
PK 2102	
BE 6701	StomaSoft Laryngectomy Tube 9/55 Non-fenestrated, 1 ea
BE 1055EF	HME Cartridge EasyFlow Heat and Moisture Exchange, pack of 30
PK 2103	
BE 6703	StomaSoft Laryngectomy Tube 10/55 Non-fenestrated, 1 ea
BE 1055EF	HME Cartridge EasyFlow Heat and Moisture Exchange, pack of 30
PK 2106	
BE 6705	StomaSoft Laryngectomy Tube 12/55 Non-fenestrated, 1 ea
BE 1055EZ	HME Cartridge ClassicFlow Heat and Moisture Exchange, pack of 30
PK 2107	
BE 6705	StomaSoft Laryngectomy Tube 12/55 Non-fenestrated, 1 ea
BE 1055EF	HME Cartridge EasyFlow Heat and Moisture Exchange, pack of 30
PK 2104	
EL 1000	Electrolarynx, 1 ea
BE 6703	StomaSoft Laryngectomy Tube 10/55 Non-fenestrated, 1 ea
BE 1055EZ	HME Cartridge ClassicFlow Heat and Moisture Exchange, pack of 30
PK 2105	
IN 1608-SO	16 Fr. 8 mm Classic Indwelling VP Sterile Option, 1 ea
IN 1610-SO	16 Fr. 10 mm Classic Indwelling VP Sterile Option, 1 ea
IN 4069	16 Fr. Plug Insert, 1 ea
BE 6703	StomaSoft Laryngectomy Tube 10/55 Non-fenestrated, 1 ea
BE 1055EZ	HME Cartridge ClassicFlow Heat and Moisture Exchange, pack of 30



PK 2105



PK 2101

Blom-Singer® Adjustable Bi-Flanged Fistula Protheses (Sterile)

Ref	Description
FM 1003*	25mm (Small) Adjustable Bi-Flanged Fistula Prosthesis, 1 ea
FM 1004*	38mm (Medium) Adjustable Bi-Flanged Fistula Prosthesis, 1 ea
FM 1005	50mm (Large) Adjustable Bi-Flanged Fistula Prosthesis without IR Catheter
FM 1010	Blom-Singer Fistula Prosthesis Insertion and Removal Catheter

*FM 1010 Insertion and Removal Catheter supplied at no charge with FM 1003 and FM 1004.

Manometer with Manometer Adaptor


A tool to teach patients how to self-regulate intratracheal air pressure. Proper airflow control may allow for better voice quality, less fatigue, and longer duration of the tracheostoma adhesive housing seal.

Ref	Description
BE 8000	Manometer Kit includes 1 manometer, 1 manometer adaptor and 1 adhesive housing
BE 8020	Replacement Manometer Adaptor includes 5 manometer adaptors and 5 adhesive housings

Blom-Singer® Insufflation Test Set

The Blom-Singer Insufflation Test Set is designed to aid the clinician in subjectively predicting the ability of the laryngectomized patient to produce tracheoesophageal voice.

Ref	Description
BE 6040	Insufflation Test Set, box of 6



The Choice for Long- or Short-Term Augmentation

RENÚ® VOICE AND RENÚ® GEL

RENÚ Voice is for long-term augmentation:

- A buffered hydrogel of synthetic, proven-biocompatible calcium hydroxylapatite (CaHA) particles 25 to 45 microns in size for consistent and easy injection.¹

RENÚ Gel is for short-term augmentation:

- A buffered hydrogel of sodium carboxymethylcellulose² (NaCMC) and glycerin.
- Resorbs slowly in vivo, providing a scaffold for local tissue infiltration.³

RENÚ® Voice: Calcium Hydroxylapatite Implant

RENÚ Voice is an injectable vocal fold implant intended for long-term augmentation.

Ref	Description
08-015-04-V04	RENÚ Voice Calcium Hydroxylapatite Implant, 1.5cc
08-015-04-VOI	RENÚ Voice Calcium Hydroxylapatite Implant, 1.5cc, CE

RENÚ® Gel: Hydrogel Implant

RENÚ Gel is an injectable vocal fold implant intended for short-term augmentation.

Ref	Description
08-015-00-V04	RENÚ Gel Hydrogel Implant, 1.5cc

RENÚ® Transoral Needle

RENÚ ENT Needle is intended to be used for the trans-oral injection of RENÚ Voice and RENÚ Gel.

Ref	Description
12-000-00-ND1	RENÚ Transoral Needle, 24 gauge needle tip, 9.8" long, 16 gauge shaft



RENÚ® VOICE



RENÚ® GEL



Needle tip

1. Data on file, InHealth Technologies®.


2. Mallur PS, et al. Safety and efficacy of carboxymethylcellulose in the treatment of glottic insufficiency. The Laryngoscope. 2012;122:322-326. <https://doi.org/10.1002/lary.21930>.

3. Hughes RGM, Morrison M. Vocal cord medialization by transcutaneous injection of calcium hydroxylapatite. Journal of Voice. 2005;19(4):674-678.

Katz Extractor®

OTO-RHINO FOREIGN BODY REMOVER

The easy way out



Katz Extractor® Oto-Rhino Foreign Body Remover

Designed in conjunction with David L. Katz, MD, MPH, the Katz Extractor Oto-Rhino Foreign Body Remover is intended for extraction of impacted foreign bodies in the nasal passage. Single use only. This device may also be used in the external ear canal, based upon clinical judgment. Covered by one or more of the following patents and other patents pending: US#5,454,817. Katz Extractor® is a registered trademark of Freudenberg Medical, LLC.

Ref	Description
FB 5000	Oto-Rhino Foreign Body Remover, box of 5




Nasal Septal Perforation Prostheses


The Blom-Singer® Nasal Septal Perforation Prosthesis is a two-piece, magnetically coupled solution for non-surgical closure of nasal septal perforations.

OVAL


Ref	Description
SP 2112	Nasal Septal Perforation Prosthesis Oval 11 x 13 x 2mm
SP 2113	Nasal Septal Perforation Prosthesis Oval 11 x 13 x 3mm
SP 2132	Nasal Septal Perforation Prosthesis Oval 13 x 17 x 2mm
SP 2133	Nasal Septal Perforation Prosthesis Oval 13 x 17 x 3mm
SP 2152	Nasal Septal Perforation Prosthesis Oval 15 x 21 x 2mm
SP 2153	Nasal Septal Perforation Prosthesis Oval 15 x 21 x 3mm




Oval 11 x 13mm



Oval 13 x 17mm




Oval 15 x 21mm




ROUND


Ref	Description
SP 1052	Nasal Septal Perforation Prosthesis Round 5 x 2mm
SP 1053	Nasal Septal Perforation Prosthesis Round 5 x 3mm
SP 1072	Nasal Septal Perforation Prosthesis Round 7 x 2mm
SP 1073	Nasal Septal Perforation Prosthesis Round 7 x 3mm
SP 1102	Nasal Septal Perforation Prosthesis Round 10 x 2mm
SP 1103	Nasal Septal Perforation Prosthesis Round 10 x 3mm




Round 5mm



Round 7mm



Round 10mm



ORDERING INFORMATION

UNITED STATES

Phone: Toll-Free: 800.477.5969

Fax: Toll-Free: 888.371.1530

E-mail: order@inhealth.com

Prescription E-mail: ihtrx@inhealth.com

Order Online: www.inhealth.com

Mail: InHealth Technologies® • 1110 Mark Ave. • Carpinteria, CA 93013 • US

INTERNATIONAL

InHealth Technologies® products are available through authorized distributors worldwide. For international distributor information, please call our customer service department at **+1.805.684.9337** or e-mail us at order@inhealth.com.

NEW PATIENT ENROLLMENT INSTRUCTIONS

Enrolling as a new patient? It's easy! Here is what to do:

1. Obtain a **prescription** from your provider.
2. Fill out the **Client/Patient Information Sheet**.
3. Send back to InHealth Technologies® with a copy of your insurance card (front and back). You can send via email: order@inhealth.com or fax: **888.371.1530**
4. Place your order! Please remember that a prescription alone will not automatically generate an order.

Need insurance or Medicare help? Contact us!

Phone: **800.477.5969 option 3**

Email: medicare@inhealth.com

While you don't have to fill anything out on the **Client/Patient Bill of Rights** form, please make sure you review it. We want you to know your rights.

All forms are available under the **Resources** tab on the InHealth Technologies® website: www.inhealth.com

FREE EDUCATIONAL RESOURCES

- Available on www.inhealth.com under the **Resources** tab.
- Visit InHealth Technologies® YouTube channel at www.cutt.ly/inhealthvideo to view our videos.

Rx: Products with the Rx notation require a prescription from a physician (US patients).

In 1978, Drs. Singer and Blom pioneered the methods of tracheoesophageal puncture and valved silicone voice prostheses that have been the international standard for voice restoration for over 40 years, allowing thousands of patients worldwide to regain their ability to speak. As a leading supplier of advanced voice restoration systems, our goal at InHealth Technologies® is to help laryngectomees regain and preserve a vital part of their lives—their voice. We are proud to continue the Blom-Singer® legacy today with clinical support, input, and leadership.

Blom-Singer®



STAY CONNECTED...

Speakers Club™

**Patients > Please visit us at www.inhealth.com/speakersclub
SMS group texting: **SPEAK to 855.909.0670****

We created this digital newsletter for you! Each month, the **Speakers Club News** will provide educational laryngectomy and voice prosthesis information.

Receive our latest updates with our SMS group texting option. **Text SPEAK to 855.909.0670** to receive news and updates for patients and caregivers.

ProNET™

**Professionals > Please visit us at www.inhealth.com/pronet
SMS group texting: **PRO to 855.909.0670****

We created this digital newsletter specifically for professionals. Each month, the **ProNET News** will deliver useful laryngectomy and voice prosthesis information for you and your practice, including material you may want to share with your patients.

Receive our latest updates with our SMS group texting option. **Text PRO to 855.909.0670** to receive announcements, new product information and educational opportunities.

INHEALTH
TECHNOLOGIES®

United States
1110 Mark Ave., Carpinteria, CA 93013
Toll-Free: 800.477.5969
Toll-Free Fax: 888.371.1530

International
Tel: +1.805.684.9337
Fax: +1.805.684.8594
E-mail: order@inhealth.com

Patents: www.inhealthpatents.com

Terms & Conditions: www.inhealth.com/terms

©2020 InHealth Technologies®—Manufactured by Freudenberg Medical, LLC (191207.00)

Blom-Singer®

The voice of experience since 1978



Order online at www.inhealth.com